VIDAS[®] TOXO Compétition (TXC)

IVD

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VIDAS[®] TOXO Competition (TXC)

IVD

VIDAS TOXO Competition is an automated qualitative test for use on the VIDAS instruments for the enzyme immunoassay detection of anti-Toxoplasma gondii total Ig in human serum or plasma (lithium heparinate or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Toxoplasma gondii, an obligate intracellular protozoan parasite, is a significant pathogen among humans. The parasite, whose primary host is the Felidae, is scattered in nature and invades all orders of mammals.

Toxoplasmosis is usually benign or asymptomatic, but can have severe consequences if it occurs in immunodeficient subjects or fetuses (2). Congenital toxoplasmosis connected with maternal contamination prior to conception is exceptional and most often occurs in immunocompromised women (7). Women who are seropositive before they become pregnant are essentially protected from transmitting the infection to their unborn child. Women who are seronegative are at risk of becoming infected during gestation (1).

The diagnosis of *Toxoplasma* infection is most commonly made by biological examination: specific immunoglobulin detection (IgM and IgG)(4, 5, 6).

VIDAS TOXO Competition is a screening test designed to assess the immunity status of patients (3). If the test is positive, seroconversion or acquired immunity should be determined using specific tests for the detection of IgG and IgM.

PRINCIPLE

The assay principle combines an inhibition-competition enzyme immunoassay method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR), serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips. All of the assay steps are performed automatically by the instrument.

The reaction medium is cycled in and out of the SPR several times.

After dilution, the sample is incubated with the SPR. Anti-Toxoplasma antibodies (IgM and IgG) present in the specimen will bind to the *Toxoplasma gondii* lysate proteins coating the interior of the SPR. Unbound components are eliminated during the wash steps.

The solid phase is then incubated with the conjugate: alkaline phosphatase labeled monoclonal anti-P30 antibodies. This conjugate will compete with the antibodies coated on the interior of the SPR by the Toxoplasma antigen. Unbound conjugate is removed by washing.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm.

The intensity of the fluorescence is inversely proportional to the concentration of immunoglobulins present in the sample. Results are automatically calculated by the instrument and expressed as an index in relation to the calibration curve stored in memory.

CONTENT OF THE KIT (60 TESTS):

| | | |
|---------------------------------------------|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 60 TXC strips | STR | Ready-to-use. |
| 60 TXC SPRs 2 x 30 | SPR | Ready-to-use. SPRs coated with <i>Toxoplasma gondii</i> lysate proteins, RH Sabin strain grown in mice (8). |
| TXC positive control 1 x 1.9 ml (liquid) | C 1 | Human serum* containing anti-toxoplasma IgG + 1 g/l sodium azide and protein stabilizers. Ready-to-use. Index: range indicated on the MLE card after the following mention : "Control C1 (+) Test Value Range". |
| TXC negative control 1 x 1.9 ml (liquid) | C 2 | Human serum* + 1 g/l sodium azide and protein stabilizers. Ready-to-use. Index: range indicated on the MLE card after the following mention : "Control C2 (-) Test Value Range". |
| Standard 1 x 1.2 ml (liquid) | S 1 | Human serum* containing anti-Toxoplasma IgG + 1 g/l sodium azide and protein stabilizers. Ready-to-use. |
| 1 MLE card | | Specifications sheet containing the factory master calibration data required to calibrate the test. |
| 1 Package insert | | |

* This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR

The interior of the SPR is coated during production with *Toxoplasma gondii* lysate proteins. Each SPR is identified by the TXC code. Only remove the required number of SPRs from the pouch and **reseal the pouch correctly after opening**.

The Reagent Strip

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the TXC strip

| Wells | Reagents |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | Sample well |
| 2 | Sample diluent: TRIS buffer (50 mmol/l) pH 7.3 - Tween + protein and chemical stabilizers + 1 g/l sodium azide (400 µl) |
| 3 | Pre-wash solution: TRIS (50 mmol/l) pH 7.3 - Tween + protein and chemical stabilizers + 1 g/l sodium azide (600 µl) |
| 4 - 5 - 7 - 8 | Wash solution: TRIS (10 mmol/l) pH 7.3 - Tween + 1 g/l sodium azide (600 µl) |
| 6 | Conjugate: Alkaline phosphatase labeled monoclonal anti-P30 antibodies (mouse) + 1 g/l sodium azide (400 µl) |
| 9 | Empty well |
| 10 | Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl). |

* IRRITANT reagent:

- **R 36:** Irritating to eyes.
- **S 26:** In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

For further information, refer to the Safety Data Sheet available on request.

MATERIAL REQUIRED BUT NOT PROVIDED

- Pipette with disposable tip calibrated to dispense 125 µl.
- Powderless, disposable latex gloves.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The optical cuvette with substrate (well 10) contains an irritant agent (6.6% diethanolamine). Refer to the risk phrase "R" and the precautions "S" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent and a solution of household bleach containing at least 0.5% sodium hypochlorite. See the Operator's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The VIDAS and mini VIDAS instruments should be regularly cleaned and decontaminated (see the Operator's Manual).

STORAGE CONDITIONS

- Store the VIDAS TXC kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS

Specimen type and collection:

Human serum or plasma (anticoagulants: lithium heparinate, EDTA). It is recommended that each laboratory checks the compatibility of collection tubes used. **The use of clearly hemolyzed, lipemic or icteric samples has not been validated.** Sera can be inactivated before testing (30 minutes at 56°C).

Specimen stability

Samples may be freshly collected or stored for 5 days at 2-8°C after which time they must be frozen at -25 ± 6 °C. Avoid successive freezing and thawing.

A study performed on frozen samples over a period of 2 months, showed that the quality of results is not affected.

INSTRUCTIONS FOR USE

For complete instructions, see the VIDAS or mini VIDAS Operator's Manual.

Master lot data entry

Before each new lot of reagents is used, specifications (or factory master calibration curve data) must be entered into the instrument (VIDAS or mini VIDAS) using the master lot entry (MLE) card (specifications sheet) included in each kit. If this operation is not performed **before initiating the tests**, the instrument will not be able to print results. The master lot data need only be entered once for each lot.

It is possible to enter data automatically using the MLE card or manually.

Calibration

Calibration, using the standard provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data has been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested in **duplicate** (see Operator's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.
2. Remove one TXC strip and one TXC SPR from the kit for each sample, control or standard to be tested. **Make sure the storage pouch has been resealed after the required SPRs have been removed.**
3. Type or select " TXC " on the instrument to enter the test code. The standard must be identified by "S1", and tested in **duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.
4. Mix the standard, controls and samples using a Vortex-type mixer.
5. Pipette **125 µl** of standard, sample or control into the sample well.
6. Insert the VIDAS SPRs and strips into the instrument. Check to make sure the color labels with the three letter assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the VIDAS Operator's Manual. All the assay steps are performed automatically by the instrument. The assay will be completed within approximately 40 minutes.
8. After the assay is completed, remove the SPRs and strips from the instrument.
9. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The instrument calculates a test value for each sample. This value is the ratio between the test's RFV (Relative Fluorescence Value) and the standard stored in memory.

Threshold and interpretation of results

| Test value | Interpretation |
|--------------|----------------|
| $i \geq 1.6$ | Negative |
| $i < 1.6$ | Positive |

If the result is negative, the patient is not immunized and there is no seroconversion.

If the result is positive, it is recommended that the same sample, or another sample taken 21 days later be tested using specific tests (for IgG and IgM detection) to determine whether the result is a seroconversion or an established immunity.

Interpretation of test results should be made taking into consideration the patient history, and any other tests performed.

The VIDAS TOXO Competition reagent is calibrated against collection sera.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS TXC kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2.

Results cannot be validated if the control values deviate from the expected values.

The expected negative control value must be included in the range indicated on the MLE card.

Note

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient history, and any other tests performed.

VIDAS TOXO Competition has not been validated for use with neonatal specimens (cord blood, etc.) nor with specimens other than human serum.

PERFORMANCE

Studies performed using VIDAS TXC gave the following results:

Sensitivity and specificity:

A multicenter study on 2021 sera including 200 sera from seroconversion, was conducted. Sera which gave a discrepant result in comparison with reference techniques were tested in Dye-Test and ISAGA. 7 sera could not be classified by reference techniques and were excluded from the performance calculation.

Results were the following:

| | | N = 2014 | |
|-----------|---|----------------------|-----|
| | | Reference techniques | |
| | | + | - |
| VIDAS TXC | + | 1101 | 6 |
| | - | 7 | 900 |

Sensitivity: 99.37 %. Range at 95% (98.68 – 99.70 %)

Specificity: 99.34%. Range at 95% (98.54 – 99.70 %)

The only serum from seroconversion which was not detected with VIDAS Toxo Competition had the following characteristics:

Dye-Test: 2 IU/ml (threshold value)

Sensitized Direct Agglutination: Negative

ISAGA: 6 (equivocal).

PRECISION

Within-run reproducibility:

3 samples were tested in a same run

| | n | RFV moyenne | CV (%) |
|-------------------|----|-------------|--------|
| negative | 30 | 2215 | 2.8 |
| weakly positive | 26 | 1700 | 2.5 |
| strongly positive | 27 | 118 | 6.2 |

Between-run reproducibility

3 samples were tested singly in 17 different runs on the same VIDAS instrument.

| | mean RFV | CV (%) |
|-----------------|----------|--------|
| negative | 2117 | 4.4 |
| weakly positive | 1706 | 3.7 |
| positive | 1044 | 4.6 |

CROSS REACTIVITY AND RELEVANT INTERFERENTS

| | | Immunologic status in relation to Toxoplasmosis according to reference techniques | |
|----------------------------|----|-----------------------------------------------------------------------------------|----------|
| | | Positive | Negative |
| HBs Ag | 11 | 9 | 2 |
| HIV | 14 | 10 | 4 |
| P. falciparum | 10 | 6 | 4 |
| CMV IgM | 8 | 7 | 1 |
| Anti-nuclear antibodies | 23 | 7 | 16 |
| EBV | 20 | 2 | 18 |
| Lyme | 4 | 2 | 2 |
| Multiparous pregnant women | 21 | 8 | 13 |
| C. trachomatis | 7 | 7 | 0 |
| Rheumatoid factors | 4 | 0 | 4 |
| Syphilis | 4 | 4 | 0 |

No interference linked to the presence of these potential interferents was revealed, except for 2 EBV-positive samples and one sample with rheumatoid factor which gave false positive results with VIDAS TOXO Competition.

PREVALENCE

T. gondii is a strict pathogen whose prevalence differs from one country to another or even one region to another.

Contamination by *T.gondii* can vary according to cultural customs and eating habits, resulting in a prevalence ranging from less than 10% in certain regions of Northern Europe to more than 90% in Africa.

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other disposable contaminated materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

LITERATURE REFERENCES

- Bertrand LECOLIER " Séroconversion de la toxoplasmose "Tempo Médical n° 422-20/03/91-13.
- FORTIER B., AJANA F., CAMUS D. "Prévention, diagnostic et suivi de la toxoplasmose congénitale". NPN, Médecine n°165-Mai 1991.
- DAVE J., BALFOUR A.H., and PERKIN J. "Evaluation of the VIDAS toxo competition assay for the detection of antibodies to *Toxoplasma gondii* in human sera". Serodiagnosis Immunotherapy Infectious Disease. 1994, 6, 41-45.
- POULETTY P., PINON J.M., GARCIA GONZALEZ M., DESONTS G., THULLIEZ P., THOANNES H., KADOUCHE J. An Anti-Human Immunoglobuline M Monoclonal Antibody for Detection of Antibodies to *Toxoplasma gondii*. Eur. J. Clin.Microbiology 1984, 3, 510 - 515.
- SANTORO F., AFCHEIN D., PIERCE R., CESBRON J. Y., OVLAQUE G., CAPRON A. Serodiagnosis of *Toxoplasma* infection using a purified parasite protein (P30), Clin. Exp.Immunol., 1985, 62, 262 - 269.
- Y. NAOT, J.S. REMINGTON. An Enzyme-Linked Immunosorbent Assay for Detection of IgM Antibodies to *Toxoplasma gondii* : Use for Diagnosis of Acute Acquired Toxoplasmosis. Jnl Inf Dis. 1980, 142, 757 - 766.
- G. DESMONTS. – Toxoplasmoses congénitales. Cinq cas de transmission à l'enfant d'une infection maternelle antérieure à la grossesse. – Press. Méd. 1990 : 1445-9.
- COUZINEAU P. and BAUFINE-DUCROCQ H. Study of the possibilities of utilization of TG 180 sarcoma of the mouse. Application to toxoplasmosis. Ann.Parasitol.Hum.Comp, 1969, 44, p.217-224.

INDEX OF SYMBOLS

| Symbol | Meaning |
|-----------------------------------------------------------------------------------|----------------------------------------------|
|  | GB : Catalogue number US : Catalog number |
|  | In Vitro Diagnostic Medical Device |
|  | Manufacturer |
|  | Temperature limitation |
|  | Use by |
|  | Batch code |
|  | Consult Instructions for Use |
|  | Contains sufficient for <n> tests |

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

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bioMérieux, the blue logo, VIDAS and SPR are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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